



# SenoBright™

## CONTRAST-ENHANCED SPECTRAL MAMMOGRAPHY

When a mammogram and ultrasound are inconclusive, SenoBright offers a straightforward, cost-effective problem-solving alternative to MRI that delivers faster, confident diagnosis for a better patient experience.

**Increase diagnosis confidence**

**HIGH SPECIFICITY**

for breast cancer diagnosis

**HIGH SENSITIVITY**

to reduce false-positives and help prevent unnecessary biopsies and surgeries

**LESS THAN 7-MINUTE EXAM**

**Improve patient experience**



**LESS WAITING**

Perform the quick exam at the same time, same location, with the same equipment



**REDUCED ANXIETY**

Release your patient immediately in the case of a negative exam

**Reduce Costs**



**LESS EXPENSIVE...**

alternative diagnostic method



**FEW FALSE POSITIVES...**

high specificity helps reduce unnecessary biopsies and surgeries



**ENHANCED MRI AVAILABILITY...**

frees up MRI time for other exams

## How It Works

**Contrast agent**

*highlights areas of unusual blood flow*

**4** simple steps in your clinical routine

**1**  
INSTALL

**2**  
INJECT

**3**  
ACQUIRE

**4**  
REVIEW

**STANDARD INTRAVENOUS IODINE INJECTION**  
performed in your radiology department

**ACQUISITION OF USUAL 4 MAMMOGRAPHIC VIEWS**  
on the same mammography system